


Hbouis Tradis 2019-03-18


→ BGD → should be BCD.



Bacterial culture after transduction



Plate onto Kanamycin to select for transductants.



1 plate (~ 5000 colonies)
scraped
resuspended in ~ 3ml

↓
Frozen in
6 x 0.5ml
aliquots

Sublibrary B



2 plates

↓
6ml media (liquid)

↓
Frozen in 12 x 0.5ml
aliquots

Sublibrary C

Plated over 10 plates
≈ 5,000 per plate

7 plates

↓ liquid.
21ml media

↓
Frozen in 42 x 0.5ml
aliquots

Sublibrary D

Sublibrary B ~ 5000 colonies
Sublibrary C ~ 10,000 colonies
Sublibrary D ~ 35,000 colonies